


NEW PROSTHETIC OPTIONS

HIGH PROSTHETIC VERSATILITY



XCN[®] IMPLANT SYSTEM



FOR NATURAL SMILES ONLY



Orthodontics and Implantology


Leone S.p.a.
Via P. a Quaracchi 50
50019 Sesto Fiorentino | Firenze | Italy
Export Dept.:
phone +39 055.3044620 | fax +39 055.304405
export@leone.it | www.leone.it



Orthodontics and Implantology

IN-06-19


**XCN® IMPLANT SYSTEM:
FOUR DIFFERENT IMPLANT MACRO-DESIGNS**



MAX STABILITY IMPLANTS

(root-form geometry, aggressive thread design)

- poor bone density
- post-extraction sockets
- immediate loading
- ridge split



CLASSIX IMPLANTS

(cylindrical geometry, atraumatic thread design)

- medium and high bone density
- limited horizontal bone availability
- crestal sinus lift



2.9 NARROW IMPLANTS

(diameter of only 2,9 mm)

- limited interdental spaces and narrow ridges
- upper lateral incisors
- lower central and lateral incisors



6.5 SHORT IMPLANT



(length of only 6,5 mm)

- limited vertical bone availability



MORSE TAPER IMPLANT-ABUTMENT CONNECTION

All XCN® implants present the key features of the Leone Implant System, especially the screwless self-locking Morse taper connection and the Platform Switching design with all the well-known advantages, outstanding mechanical strength and stability, no micro-gaps and micro movements at the implant-abutment interface, preservation of crestal bone over time and prosthetic simplicity and safety.





NEW PACKAGINGS



NEW IMPLANT CARRIER

Pre-mounted, colour coded, built-in torque stop, 2 depth marks, easy removal



**PROSTHETIC SIMPLICITY:
ONLY 2 IMPLANT-ABUTMENT CONNECTION**

NEW KITS

Surgical, Osteotomy, Depth Stops, Organizer



NEW DRILLS

With DLC coating: high resistance to wear and corrosion, excellent visibility of drilling depth

